

# Short Manual for the MRI Measurement with multi-parameter mapping protocol for Siemens 3T MRI


## 1. Standard Operation Procedure

### 1.1 Explanation of the Required Sequences and How-to-Measure

The MR protocol is composed of the following sequences (parts A-D).

#### A. Localizer

Please place the `localizer` in order to cover the whole brain and also the cervical cord until C2/C3 level for the multi-parameter mapping (MPM) protocol (see below).



#### B. RF-Transmit-Map (B1 Mapping)

The sequence `rf_map` records a map of the RF transmit field. This sequence can be acquired directly after the localizer before running the MPM protocols, applying the same planning, i. e. whole brain coverage.

#### C. RF-Sensitivity-Maps

For obtaining the RF sensitivity field, we apply two separate sequences right before each multi-parametric mapping sequence (see also next section). One is `RF_sens_head` and the other is `RF_sens_body`, whereas for the latter the Body Coil has to be activated in the System Card as indicated below:



**Attention:** Although the automatic coil selection should be deactivated, the process of copying the protocol parameters manually or via copy reference (“Center of slice groups & sat. regions”) might also copy the coil selection will be copied. Thus you have to **check and change the coil selection** accordingly whenever you manually copy or set a copy reference from RF\_sens\_body to another sequence and vice versa!

## D. Multi-parametric mapping (MPM) sequences

This protocol is considered as a research protocol in the frame of the NISCI trial and will be processed to obtain the quantitative MR parameters such as magnetization transfer (MT), and relaxation parameters R1 and R2\*. These parameters are sensitive to the myelination changes and iron level in the tissue.

The protocol comprises of following three multi echo 3D FLASH sequences:

- T1-weighted sequence
- PD-weighted sequence
- MT-weighted sequence

The order is not important within the MPM protocol, but for each scan the RF\_sens\_head and RF\_sens\_body sequences have to be run before each of the three measurements.

Please make use of copy references in order to have the same acquisition volume for the T1-, PD- and MT-weighted sequences, while checking for the correct coil usage (see above).

## 2. Listing of all Sequences within the Study Protocol

1. Head Localizer
  2. RF\_map
  3. RF\_sens\_head
  4. RF\_sens\_body
- } These two sequences will be repeated before each 3D FLASH sequence of the MPMs.
5. MPM protocol with iPAT3/iPAT2 consisting of three 3D FLASH sequences, i. e. MT-, PD- and T1-weighted sequences covering the whole brain